


# Textured Polyester & Nylon Multifilament Yarns

## For Embolic Coil Device Applications



*The textured yarns are heat-set and scoured to produce an open yarn bundle.  
The softly crimped filaments blossom to aid in the occlusion of blood vessels.*

Item	Titer <sup>1</sup> (Denier)	Filaments	Tenacity <sup>2</sup> (g/d)	Elongation (%)	Net Wt. (oz)	Package Type
40/34-PA66-FTT-SS	40	34	4.2	28	~12	Cone
44/27-PET-5540-FTT-SS	40	27	4.2	32	~16	Cone

1. Titer Conversion:  $Dtex = Denier \times 1.1111$
2. Tenacity Conversion:  $cN/dtex = g/d \times 0.8826$

**DISCLAIMER:** The information supplied in this document is for guidance only and should not be construed as a warranty. All implied warranties are expressly disclaimed, including without limitation any warranty of merchantability and fitness for use. All users of the material are responsible for assuring that it is suitable for their needs, environment and end use. All data is subject to change as TDA deems appropriate.